

Verification of BioFire® FilmArray Pneumonia Panel *plus* Detections by Alternate Methods

P0569
Abstract #2752


D. Judd¹, B. Graham¹, D. Sadrija¹, O. Cham¹, C. Li¹, J. Stone¹, K. Broadbent¹, C. Graue¹, M. Jones¹, J. Cloud¹, and M. Rogatcheva¹

¹BioFire Diagnostics, LLC, Salt Lake City, UT, USA.

Background

The BioFire FilmArray Pneumonia Panel *plus* is intended to identify pneumonia causing agents in unprocessed sputum (including endotracheal aspirate (ETA)) and bronchoalveolar lavage (BAL) specimens. The panel detects bacteria, viruses, and select antimicrobial resistance markers. Clinical performance was evaluated in a multi-center study by comparing FilmArray detections to the results of other methods used to detect pathogenic organisms.

Any discordant detections were investigated to determine the root cause of false positive (FP) and false negative (FN) discrepancies using additional molecular methods.



BioFire® FilmArray® Pneumonia Panel *plus* Investigational Use Only
34 Targets

Bacteria <i>Acinetobacter calcoaceticus-baumannii</i> complex <i>Serratia marcescens</i> <i>Proteus</i> spp. <i>Klebsiella pneumoniae</i> group <i>Enterobacter aerogenes</i> <i>Enterobacter cloacae</i> complex <i>Escherichia coli</i> <i>Haemophilus influenzae</i> <i>Moraxella catarrhalis</i> <i>Pseudomonas aeruginosa</i> <i>Staphylococcus aureus</i> <i>Streptococcus pneumoniae</i> <i>Klebsiella oxytoca</i> <i>Streptococcus pyogenes</i> <i>Streptococcus agalactiae</i>	Atypical Bacteria <i>Legionella pneumophila</i> <i>Mycoplasma pneumoniae</i> <i>Chlamydia pneumoniae</i> Viruses Influenza A Influenza B Respiratory Syncytial Virus Human Rhinovirus/Enterovirus Human Metapneumovirus Parainfluenza Virus Adenovirus Coronavirus Middle East Respiratory Syndrome Coronavirus	Antimicrobial Resistance Genes <i>mecA/C</i> and <i>MREJ</i> KPC NDM OXA-48-like CTX-M VIM IMP
---	---	--

Sample Type: Sputum, Endotracheal aspirate, Bronchoalveolar Lavage, and mini-BAL

BioFire Pneumonia Panel binning algorithm used to report semi-quantitative results

Assay Result	Reported Result and Bin
Negative or <10 ^{3.5} copies/mL	Not Detected
Positive and ≥10 ^{3.5} – <10 ^{4.5} copies/mL	Detected at 10 ⁴ copies/mL
Positive and ≥10 ^{4.5} – <10 ^{5.5} copies/mL	Detected at 10 ⁵ copies/mL
Positive and ≥10 ^{5.5} – <10 ^{6.5} copies/mL	Detected at 10 ⁶ copies/mL
Positive and ≥10 ^{6.5} copies/mL	Detected at ≥10 ⁷ copies/mL

Materials and methods

All comparator testing was performed in a blinded manner by personnel with no knowledge of BioFire Pneumonia Panel results

Comparator Methods:
Quantitative reference culture (qRefCx) is a standard culture method that relies on plate enumeration to obtain a quantity. qRefCx was the comparator method used for semi-quantitative bacterial detections and was performed at a central reference laboratory. Results quantified at or above 10^{3.5} CFU/mL were considered positive while results below the cutoff were considered negative.

Quantitative molecular method (qMol) is PCR followed by Next-Generation Sequencing (NGS). qMol was the comparator method used for AMR detections and was performed at the BioFire Laboratory.

Conventional PCR and Sanger sequencing was used in this study and utilizes 2 specific molecular assays for all atypical bacteria and viral analytes. Conventional PCR followed by Sanger sequencing was the comparator method used to confirm viral and atypical bacterial detections and was performed at the BioFire Laboratory.

Clinical samples with discordant results between FilmArray and comparator testing were selected for investigation and DNA from residual specimens was extracted and tested using molecular assays targeting genes or gene regions different from the FilmArray targets. Bi-directional sequencing of the amplified region(s) was used to confirm the identity of organisms.

Discrepancy Methods:
Standard of Care (SOC) testing is a culture method ordered by the physician and performed at the medical clinic when a patient presents with symptoms of a lower respiratory tract infection. SOC results were used as an initial method for resolving discordant detections between FilmArray and the comparator method.


Other molecular methods refers to additional PCR based molecular assays followed by Sanger sequencing. These additional assays were designed at BioFire Diagnostics and performed in the BioFire Laboratory.

FilmArray Pneumonia Panel retests using residual sample were also used to resolve discrepancies. The retests were performed at the BioFire Laboratory and were mainly used when investigating false negative results.

Bacteria

Confirmation of BioFire FilmArray Detections by Different Methods


Reference Method: qRefCx ≥10^{3.5} CFU/mL



Viruses and Atypical Bacteria

Confirmation of BioFire FilmArray Detections by Different Methods


Reference Method: Conventional PCR and Sanger Sequencing



Antimicrobial Resistance (AMR) Genes



Confirmation of BioFire FilmArray Detections by Different Methods

Reference Method: qMol




BioFire FilmArray Unbinned Values Reported

Culture Positive vs Culture Negative



BioFire FilmArray Unbinned Values Compared To qMol and qRefCx in Sputum



Histogram distribution of BioFire FilmArray unbinned values observed in samples that were qRefCx positive compared to samples that were qRefCx negative. As expected, BioFire FilmArray is more sensitive to low level analytes than qRefCx. Unexpectedly, BioFire FilmArray also detected high titer organisms that qRefCx missed.

Certain organisms are more likely to be identified by qRefCx than others. For example, *P. aeruginosa* is identified by qRefCx the majority of the time. On the other hand, the BioFire Pneumonia Panel identified *H. influenzae* more often than qRefCx.

Bacteria False Negatives

Cause of False Negative Results

Analyte	BioFire FilmArray Below Quant (low level)	Missed Detection (low level)	Misidentified
<i>E. cloacae</i> complex	1	1	0
<i>E. coli</i>	1	0	0
<i>K. pneumoniae</i> group	1	1	0
<i>S. marcescens</i>	1	0	0
<i>S. aureus</i>	1	1	0
<i>E. aerogenes</i>	0	1	1
<i>H. influenzae</i>	0	1	1
<i>P. aeruginosa</i>	1	0	2
ACB complex	0	0	1

Reference Lab Identification	16S Sequencing Results and ID Confidence	In House VITEK Results and ID Confidence	High Resolution Molecular Method Sequence Results and ID Confidence
ACB complex	99% <i>Pseudomonas</i> spp.	96% <i>P. fluorescens</i>	N/A
<i>E. aerogenes</i>	100% <i>Hafnia</i> spp.	99% <i>Hafnia alvei</i>	98% <i>Hafnia paralvei</i>
<i>H. influenzae</i>	100% <i>H. influenzae</i> or <i>Haemophilus</i> spp.	Unidentified	97% <i>H. haemolyticus</i>
<i>P. aeruginosa</i>	100% <i>P. aeruginosa</i> or <i>Pseudomonas</i> spp.	Low Discrimination-- <i>Pseudomonas</i> spp.	94% <i>P. denitrificans</i>
<i>P. aeruginosa</i>	100% <i>P. aeruginosa</i> or <i>Pseudomonas</i> spp.	Low Discrimination-- <i>Acinetobacter</i> or <i>Pseudomonas</i> spp.	99% <i>P. fluorescens</i>

These isolates were investigated and determined to be misidentified by the reference laboratory. Resolved identifications are highlighted in green.

Nature of Discrepancies

Discrepancy investigation demonstrated that the majority of discordant results were a consequence of analyte levels at or below the LoD of both the BioFire FilmArray and comparator assays (which had equivalent sensitivity) and thus resulted in spurious detections.


- Nearly all discrepancies observed in the clinical trial were sequence confirmed to be present in the sample.
- When multiple replicates per specimen were tested inconsistent detections were observed for both FilmArray and comparator assays.

Viral and Atypical Bacteria False Negatives

Analyte	Missed Detection (low level)	Unresolved
<i>L. pneumophila</i>	1	0
Adenovirus	4	0
Coronavirus	6	1
hMPV	1	0
HRV/EV	2	0
Influenzae B	1	0
Parainfluenza Virus	3	0
<i>M. pneumoniae</i>	0	1

Coronavirus detections in sputum for BioFire FilmArray and comparator assays are shown in the accompanying graph. The red line depicts the level at which Coronavirus positive samples at LoD are detected. The Coronavirus false positive results are all detected later than what would be expected at LoD thus suggesting that these samples are Coronavirus positive but at a concentration below LoD. As expected, false negative results are below the LoD of the comparator assays.

Detected Coronavirus Cp Distribution



Discrepancy Investigation Summary

	Bacteria Discrepancy Analysis		BAL	Sputum
	FN	FP		
FN	Total Discrepancy	3	25	28
	FilmArray below the LoD bin	1	1	2
	FilmArray negative; evidence by qMol	1	2	3
	Resolved-other (e.g. isolate misidentification, SOC)	1	3	4
	Not Resolved	0	1	1
FP	Total Discrepancy	328	547	875
	qRefCx enumerated below 10 ^{3.5} CFU/mL	105	115	220
	qRefCx negative; evidence by qMol	218	421	639
	Resolved-other (e.g. SOC)	3	9	12
	Not Resolved	2	2	4

15/16 (93.8%) of bacterial false negatives were resolved as low level or misidentification

871/875 (99.5%) of bacterial false positive results were resolved as low level true positive

	Atypical Bacteria Discrepancy Analysis		BAL	Sputum
	FN	FP		
FN	Total Discrepancy	0	2	6
	Evidence by independent molecular method	0	2	4
	Resolved-other	0	0	0
	Not Resolved	0	0	0
	Total Discrepancy	2	0	2
FP	Total Discrepancy	0	2	10
	Evidence by independent molecular method	0	2	6
	Resolved-other	0	0	0
	Not Resolved	0	0	0
	Total Discrepancy	2	2	4

13/15 (86.7%) of atypical bacteria and AMR false negatives were resolved as low level true positive

	Virus Discrepancy Analysis		BAL	Sputum
	FN	FP		
FN	Total Discrepancy	8	10	18
	Evidence by independent molecular method	5	10	15
	Resolved-other	0	0	0
	Not Resolved	3	0	3
	Total Discrepancy	31	33	64
FP	Total Discrepancy	19	21	40
	Evidence by independent molecular method	19	21	40
	Resolved-other (e.g. SOC, FilmArray retest)	2	3	5
	Not Resolved	0	0	0
	Total Discrepancy	21	24	45

15/18 (83.3%) of viral false negatives contained the analyte of interest at or below LoD

45/64 (70.3%) of viral false positive results were resolved as low level true positive

	Antimicrobial Resistance Gene Discrepancy Analysis		BAL	Sputum
	FN	FP		
FN	Total Discrepancy	7	6	13
	Evidence by independent molecular method	7	5	12
	Resolved-other	0	0	0
	Not Resolved	0	1	1
	Total Discrepancy	8	15	23
FP	Total Discrepancy	6	10	16
	Evidence by independent molecular method	6	10	16
	Resolved-other (e.g. qRefCx, AST, SOC, AST)	0	0	0
	Not Resolved	2	2	4
	Total Discrepancy	2	2	4

19/25 (76.0%) of atypical bacteria and AMR false positive results were resolved as low level true positive

Antimicrobial Resistance Gene Reporting Algorithm

AMR Gene Result	Applicable Bacteria
<i>mecA/C</i> and <i>MREJ</i>	<i>Staphylococcus aureus</i>
CTX-M	<i>Acinetobacter calcoaceticus-baumannii</i> complex
IMP	<i>Enterobacter aerogenes</i>
KPC	<i>Enterobacter cloacae</i> complex
NDM1	<i>Escherichia coli</i>
OXA-48-like	<i>Klebsiella pneumoniae</i>
VIM	<i>Pseudomonas aeruginosa</i> <i>Proteus</i> spp. <i>Serratia marcescens</i>

AMR genes are only reported if a host bacteria is also detected above the quantification cutoff

In depth summary of MRSA detection in the BioFire FilmArray Pneumonia Panel is presented on poster P0561 by B. Graham et al.

Conclusions

The bacterial assays in the BioFire Pneumonia Panel are highly specific and more sensitive than traditional culture methods. Across all bacterial assays for both sample types, the specificity of the panel is 96.5% compared to culture. Discrepancy investigation demonstrated that nearly all false positive detections (871/875) were correctly reported by the BioFire Pneumonia Panel. If sequence confirmed false positive detections are taken into account the specificity for all bacterial assays would be 99.9%. BioFire Pneumonia Panel exhibits higher specificity than culture due to the fundamental differences between molecular methods and culture (detects viable organisms).

The BioFire Pneumonia Panel is a sensitive and specific in vitro diagnostic device for the detection of atypical bacteria, viruses, and antimicrobial resistance genes from BAL and sputum specimens compared to molecular methods. The combined specificity across both sample types for these assays is 99.6%. Similar to what is observed with the bacterial assays, the majority of false positive detections (64/89) were correctly reported by the BioFire Pneumonia Panel. If these sequence confirmed false positive detections are taken into account the combined specificity for atypical bacteria, viral, and AMR assays would be 99.9%.

Taking discrepancy resolution into consideration, the combined specificity for all assays in the BioFire Pneumonia Panel would reach 99.9%, thus providing high confidence in detections of targeted bacteria and viruses present in lower respiratory specimens.

Data presented are from assays that have not been cleared or approved by US FDA or other regulatory agencies for diagnostic use

Contact Information:
David Judd
BioFire Diagnostics, LLC
david.judd@BioFireDx.com

