



BIOFIRE® SPOTFIRE® PRE-INSTALLATION CHECKLIST

Specifications

Power Requirements

Configuration	Voltage	Frequency	AC Power at 100 VAC		AC Power at 240 VAC	
			Active	Apparent	Active	Apparent
1 Module	100 – 240 VAC	50 – 60 Hz	117 W	121 VA	101 W	148 VA
2 Modules			188 W	192 VA	173 W	216 VA
3 Modules			259 W	264 VA	237 W	277 VA
4 Modules			328 W	333 VA	304 W	341 VA

Space: 8.52 in (W) x 12.77 in (D) (21.64 cm x 32.43 cm). Height is dependent on the number of modules as seen below:



Dimensions and Weight

Dimensions	Control Station only (W x D x H): 8.52 x 12.77 x 3.82 in (21.64 x 32.43 x 9.70 cm) Modules only (H): 4.13 in (10.49 cm) Max height (4 modules): 21.0 in (54 cm)
Weight	Maximum total weight: Approximately 74 lbs (33.6 kg) Control Station: 18 lbs (8.2 kg) Modules: 14 lbs (6.4 kg) each

Module #	Total Instrument Weight
1	32 lbs (14.5 kg)
2	46 lbs (20.9 kg)
3	60 lbs (27.2 kg)
4	74 lbs (33.6 kg)

Operation Specifications

Operations Specification	15° C – 30° C @ 15 – 80% relative humidity (non-condensing) -50 ft – 10,000 ft (-16 m – 163,048 m) operating altitude Indoor use only
---------------------------------	---

WiFi Connectivity (if desired)

- Network name & password known

Ethernet Connectivity (if desired)

- Port is live/hot
- Instrument will be located within 5 feet from a live ethernet port

Printer (if desired)

- Can be ordered separately (directly from bioMérieux)
- Connection type (USB, Ethernet, or WiFi)

Uninterruptable Power Supply (UPS) is recommended but not required

Testing Supplies

BIOFIRE® SPOTFIRE® Respiratory Panel Reagent Kit (424461)

BIOFIRE® SPOTFIRE® Respiratory Panel Mini Reagent Kit (424589)

SPOTFIRE® Quality Control Material

- SPOTFIRE® RSP Positive & Negative Controls (M425) **or to order boxes separately:**
- SPOTFIRE® RSP Positive Control (M42638)
- SPOTFIRE® RSP Negative Control (M42738)

Sample Collection Materials

- Nasopharyngeal swabs
- 1 – 3 mL transport media

Additional Testing Supplies

Bleach (10% solution) or Pre-diluted Bleach Wipes

Distilled Water

Gloves

Biohazard Bin

Biohazard Sharps Container

Regular Waste Bin

Dry Cloth or Paper Towels